

CO2 SENSOR



- Do not submerge in liquid.
- Clean with 70% isopropyl alcohol or 10% bleach solutions, disinfectants such as Steris Coverage Spray HB, ammonia, or mild soap and water.
- Rinse with water using a clean, wet, soft cloth.
- Check the sensor windows are clean and dry before using.

Sterilization

⚠ Caution

- Do not sterilize the sensor by any means.

📄 Note

- The information provided for cleaning and sterilization is offered as a guide. It may be adjusted to the protocols of each service, following the guidelines above.
- Using products and procedures that do not comply with the guidelines included in this document will accelerate the normal deterioration of the accessory, which may be damaged or malfunction.

CO2 ADAPTER



Clean with warm soap solution.

Then, submerge in a liquid disinfectant, such as:

- 70% isopropyl alcohol
- 10% sodium hypochlorite in water
- 2.4% glutaraldehyde, such as Cidex Plus
- Ammonia

Rinse with plenty of sterile water and let dry.

Sterilization

- Autoclave - only adapters for adult patients
- Ethylene oxide (ETO) - 1.5-hour exposure
- Cidex Plus solution - 10-hour soak
- Perasafe solution - 10-hour soak



Warning

Before using the adapter, check that windows are dry and free from residue, and that it has not been damaged during handling or cleaning/disinfecting.