



**neumovent ts**

Comprehensive monitoring  
High-flow oxygen therapy  
Spontaneous breathing test  
New ventilation modes  
Neonatal category



Mechanical ventilator that provides  
invasive and non-invasive ventilation  
for adult and pediatric patients,  
and for neonatal patients (Optional).

**tecme**  
committed to life

## Comprehensive monitoring

It is possible to monitor the most important parameters, from the beginning of the therapy, with P/V Flex and the Stress Index option, as well as the intermediate phases towards weaning, with P0.1 and Pimax, among other tests and variables.

## Optional modules:

### High-flow oxygen therapy

The addition of this method expands the use of the ventilator to the stages that are prior and after the mechanical ventilation itself, by providing non-invasive assistance to oxygenation.

### Spontaneous breathing test

This new tool allows a spontaneous breathing test to be performed before weaning, monitored in real time and completely safe with dedicated alarms.

### New ventilation modes

Pressure Regulated Volume Control (PRVC), Synchronized Intermittent Mechanical Ventilation with Pressure Regulated Volume Control SIMV (PRVC) and Volume Support Ventilation (VSV).

### Neonatal category

It incorporates ventilation methods that would allow patient ventilation from 5 ml tidal volume.



- Ventilation for adult and pediatric patients.
- Software upgradable for:
  - Neonatal ventilation.
  - High-flow oxygen therapy.
  - Spontaneous Breathing Test (SBT) and Stress Index.
  - PRVC, SIMV (PRVC) and VSV ventilation modes.
- Available for invasive and non-invasive ventilation.
- Automatic leak compensation.
- Comprehensive respiratory mechanics module.
- Intra-hospital transportation.
- 72 hours of trend storage.
- Low cost maintenance.
- Built-in battery with capacity higher than 2.5 hrs.
- Numerical screen.
- 2-year warranty.

50 years of innovation and development in mechanical ventilators.

[www.tecmeglobal.com](http://www.tecmeglobal.com)

#### Córdoba

Av. Circunvalación (N-O) - Agustín Tosco 3040  
X5008HJY - Córdoba, Argentina  
Tel.: +54 351 4144600 - Fax: +54 351 4144605  
e-mail: [info@tecme.com.ar](mailto:info@tecme.com.ar)

#### Buenos Aires

Trole 557 - 2º Piso, Nueva Pompeya, C.A.B.A.,  
C1437DKK - Buenos Aires, Argentina  
Tel.: +54 11 1561387326  
e-mail: [info@tecme.com.ar](mailto:info@tecme.com.ar)

#### Norcross

2825 Pacific Drive, Suite B  
Norcross, GA 30071 - USA  
P: +1 770 409 9172 - F: +1 770 729 8176  
e-mail: [info@tecme-usa.com](mailto:info@tecme-usa.com)