

**vh 2100**



## Humidifier

POWER   
Designed to deliver a comfortable level of humidity.

The VH-2100 Humidifier is part of a complete humidification system designed to deliver a comfortable level of humidity. Adjustable between nine temperature and humidification levels.

This arrangement simplifies setup, minimizing the connections required.

## Product specifications

	VH-2100-115	VH-2100-230
Dimensions	W135 x D170 x H156 mm	
Weight	1.5 Kg (w/o Chamber fitted)	
<b>Electrical Power</b>		
Supply Voltage	115 ± 12V	230 ± 23V
Supply Current	1.3 A Max.	0.7 A Max.
Supply Frequency	50/60 Hz	
Heater plate capacity	150W max	
<b>Security protection</b>		
Overcurrent protection	115 V/3 A (fast burn type)	230 V/1.5 A (fast burn type)
Overheat protection	115±3° C, thermostat.	
<b>Temperature setting</b>		
	Knob: Scale 1 to scale 9	

50 years of innovation and development in mechanical respirators.  
[www.tecmeglobal.com](http://www.tecmeglobal.com) / [info@tecmeglobal.com](mailto:info@tecmeglobal.com)

### Córdoba

Av. Circunvalación (N-0) - Agustín Tosco 3040  
X5008HJY - Córdoba, Argentina  
P: +54 351 4144600 - F: +54 351 4144605

### Buenos Aires

Trole 557 - 2º Piso, Nueva Pompeya, C.A.B.A.,  
C1437DKK - Buenos Aires, Argentina  
P: +54 11 1561387326

### Atlanta

2825 Pacific Drive. Suite B  
Norcross, GA 30071 - USA  
P: +1 770 409 9172 - F: +1 770 729 8176